

MEDICAL DEVICE POWDER-FREE SUPREME PURPLE NITRILE GLOVE LUNA



0161



EN 420:2003+A1:2009
EN 455

EN ISO 374-1:2016/A1:2018 / Type C



K

EN ISO 374-5:2016



VIRUS



- > Nitrile glove in purple color
- > Brand: Luna
- > Interior: Powder free
- > PPE Category III
- > Ambidextrous
- > Disposable
- > Sizes S, M, L and XL
- > Surface: Microtexturized fingers
- > Free of Thiuranes, Thiureas, Thiazoles, and Mercaptobenzothiazoles
- > Cuff: Rolled
- > Length: 240 mm (approx.)
- > Thickness: 0,07 ± 0,02 mm
- > Latex free
- > AQL 1,5
- > Great softness and flexibility
- > Weight (±0,3g): Size S: 3,3g / M: 3,5g / L: 3,8g / XL: 4,1g

USES AND APPLICATIONS



Cleaning



Industry



Food use



Automotion



Health sector



REGULATIONS

Regulation (EU) 2016/425, according to Personal Protective Equipment.

- ✓ EPI Category III.
- ✓ EN 420:2003+A1:2009, protective gloves, general requirements.
- ✓ EN ISO 374-1:2016/A1:2018 TYPE C, chemical protective gloves: K – Sodium Hydroxide 40%.
- ✓ EN ISO 374-5:2016, protection against microorganisms (virus, bacteria and fungi).

Regulation (EU) 2017/745, according to Medical Devices. Class I. EN 455.

- ✓ EN 455-1:2020. Medical gloves for single use - Part 1: Requirements and tests to determine the absence of holes.
- ✓ EN 455-2:2015. Medical gloves for single use - Part 2: Requirements and tests for the determination of physical properties.
- ✓ EN 455-3:2015. Medical gloves for single use - Part 3: Requirements and tests for biological evaluation.
- ✓ EN 455-4:2009. Medical gloves for single use - Part 4: Requirements and testing for shelf life determination.

IN ADDITION COMPLIES:

- ✓ UNE-EN ISO 10993-10:2011, tests for irritation and skin sensitization .
- ✓ ASTM F1671/F1671M-13 / EN 16523--:2015 (EN ISO 374-1:2016), Standard Test Method For Resistance Of Materials Used In Protective Clothing To Penetration By Blood-Borne Pathogens Using Phi-X174 Bacteriophage Penetration As A Test System.



RECOMMENDATIONS

- ✓ It is recommended to check and verify the glove before use and periodically to ensure that no damage is present, as it can significantly affect the expected protection.
- ✓ Replace immediately if any deterioration is observed in the glove.
- ✓ For conservation and storage, it must be away from sources of heat and light bulbs. Store in a cool dry place.
- ✓ It is recommended to save the inner box for possible consultation.



LOGISTICS



PALLET

- Contains:**
- 96 Cartons
 - 960 Inner boxes

LAYERS/PALLET:8
CARTONS/LAYER: 12


Width: 800 mm
Length: 1200 mm
Height: 1860 mm



CARTON

- Contains:**
- 10 Inner boxes

Width: 230 mm
Length: 315 mm
Height: 210 mm



INNER BOX

- Contains:**
- 100 Gloves

Width: 110 mm
Length: 200 mm
Height: 60 mm

SIZE	REFERENCE	INNER BOX	CARTON
S	GNV1	8436023893176	8436023893213
M	GNV2	8436023893183	8436023893220
L	GNV3	8436023893190	8436023893237
XL	GNV4	8436023893244	8436023893251

REV 12

omare.com

P.I. El Saladar II. Calle Pontevedra S/N 30564 Lorquí (Murcia). Tel: 968 68 79 50

